

INSTRUCTION MANUAL

BLACKBODY RADIATION SOURCE

Model CS1250-100

and

TEMPERATURE CONTROLLER

Model 2500E

ELECTRO OPTICAL INDUSTRIES, INC.

859 Ward Drive

Santa Barbara, CA 93111

Phone: (805) 964-6701 Fax: (805) 967-8590

OCTOBER 2004

TABLE OF CONTENTS

| <u>Section</u> | <u>Page</u> |
|---|-------------|
| INTRODUCTION | 1 |
| 1.1 DESCRIPTION | 1 |
| 1.2 LIST OF ITEMS FURNISHED..... | 1 |
| 1.3 MODEL NUMBER SPECIFICATION | 1 |
| SPECIFICATIONS | 2 |
| INSTALLATION | 3 |
| 3.1 INITIAL INSPECTION | 3 |
| 3.2 INSTALLATION | 3 |
| 3.3 REPACKING FOR SHIPMENT OR STORAGE | 3 |
| 3.4 STORAGE..... | 4 |
| OPERATION..... | 5 |
| 4.1 MODEL 2500 series CONTROLS..... | 5 |
| 4.1.1 Front Panel Controls and Indicators | 5 |
| 4.1.2 Keypad..... | 5 |
| 4.1.3 Error Message Display | 6 |
| 4.2 BASIC OPERATING PROCEDURE IN LOCAL CONTROL | 6 |
| 4.2.1 Power ON | 6 |
| 4.2.2 Setting Temperature | 6 |
| 4.2.3 Cool Down Procedure | 7 |
| 4.3 Remote Control (Option /R RS232 Interface)..... | 7 |
| THEORY OF OPERATION..... | 12 |
| 5.1 GENERAL | 12 |
| 5.2 TEMPERATURE MEASUREMENT AND CONTROL (A2)..... | 12 |
| MAINTENANCE | 13 |
| 6.1 PREVENTIVE MAINTENANCE | 13 |
| 6.2 TROUBLESHOOTING | 13 |
| CALIBRATION..... | 14 |
| 7.1 INTRODUCTION..... | 14 |
| 7.2 RECOMMENDED EQUIPMENT FOR CALIBRATION | 14 |
| 7.3 USE OF CALIBRATION DATA SHEETS (Refer to Figure 7.1)..... | 14 |
| 7.4 CALIBRATION SET-UP PROCEDURE | 17 |
| 7.5 CALIBRATION PROCEDURE | 17 |

| | |
|-------------------------------|----|
| REPLACEABLE PARTS | 20 |
| 8.1 ORDERING INFORMATION..... | 20 |

WARRANTY

APPENDIX
Certification

ILLUSTRATIONS

| <u>Figure</u> | <u>Title</u> | <u>Page</u> |
|---------------|------------------------------------|-------------|
| Figure 7.1 | Calibration Data Sheet | 16 |
| Figure 7.2 | Calibration Equipment Set Up | 19 |
| Table 7 | Type S Thermocouples | Appendix |

SECTION 1

INTRODUCTION

1.1 DESCRIPTION

The EOI source you are now using is one of a complete family of sources manufactured by Electro Optical Industries.

The EOI Model 2500 series Temperature Controller is included with the source. The temperature setting is changed via the front panel keypad, or by a remote message sent over the optional serial interface. The source temperature and set point are displayed on the front panel.

The microprocessor-based controller uses carefully selected components to provide extremely stable and repeatable control of the source. The controller is designed for minimum susceptibility to line and environmental changes.

1.2 LIST OF ITEMS FURNISHED

- (1) Blackbody Radiation Source Assembly
- (2) Model 2500 series Temperature Controller
- (3) Cables (2)
- (4) Instruction Manual

1.3 MODEL NUMBER SPECIFICATION

2500 / RS232
 └───┘ └───┘
 Output Power (E)

Add /RS232 to model number for optional RS232 serial interface.

SECTION 2

SPECIFICATIONS

EOI Blackbody Radiation Source Model CS1250-100

EOI Temperature Controller Model 2500E

| | |
|---|---|
| Temperature range | 50 °C to 1250°C |
| Stability | ± 0.25V |
| Sensing element | Platinum thermocouple (Type S) |
| Type of control | Zero crossing solid state relay, PWM control, AC Power. |
| Cavity diameter | 1.00 inch, 15 degree recessed cone |
| Cavity emissivity | 0.99 (± 0.01) |
| Warm up time to 1250°C | 50 minutes or less |
| Maximum input power to source (warm up) | 500 watts |
| Power for control at 1250°C (source) | 300 watts |
| Housing temperature | Less than 50°C at 1250°C |
| Power requirement | 100/120/200/240 VAC, 50/60 Hz, 1 phase (voltage selection by internal jumpers) |

SECTION 3

INSTALLATION

3.1 INITIAL INSPECTION

This instrument was carefully inspected before shipment.

Shipping cartons should be inspected for damage before and after unpacking. Inspect instruments for physical damage incurred in transit or while in storage. If damages have occurred during shipping, notify the carrier and Electro Optical Industries, Inc. at once. Inspect package for items and accessories as listed in Section 1.2.

3.2 INSTALLATION

- (1) Install source and temperature controller in suitable location.
- (2) Connect cable from source assembly to the circular connector on rear of temperature controller.
- (3) Connect thermocouple cable to the connector on the rear of the source assembly and to the connector on the rear of the temperature controller.
- (4) Connect temperature controller power cord to correct source of AC power.

CAUTION

Connect only to AC line power as shown on tag at the rear of instrument. Other voltages require wiring changes as shown on interconnect drawing.

3.3 REPACKING FOR SHIPMENT OR STORAGE

Using the original cartons to ship or store these instruments is not recommended. Use new cartons and follow the procedure below.

- (1) Wrap instruments in heavy plastic or other waterproof material before placing in inner container.
- (2) Fill space around instruments with packing material. Protect front panel to prevent damage to controls.
- (3) Place in outer container with packing between the containers. Seal with heavy tape or metal bands.
- (4) Mark containers "DELICATE INSTRUMENTS, FRAGILE".
- (5) DO NOT ship via parcel post.

3.4 STORAGE

If instruments are to be stored, follow packaging instructions in Section 3.3. Locate them in an approved area for storing electronic instruments.

Storage temperature range: -40°C to +50°C.

SECTION 4

OPERATION

4.1 MODEL 2500 series CONTROLS

4.1.1 Front Panel Controls and Indicators



POWER switch

Turns line power ON [1] or OFF [0]. Front panel displays illuminate when power is on.

EC Display

4 character red display indicates the current source temperature or error codes.

SET Display

4 character green display above keypad indicates the current temperature setting controlled by the  and  keys.

4.1.2 Keypad



Keypad



MENU ACCESS KEY

Used to enter the menu system, index to the next menu and enter the security level menu.



RAISE KEY

Used to increase selected parameter value and set point temperature. When held, value increases rapidly after a short delay.



LOWER KEY

Used to decrease selected parameter value and set point temperature. When held, value increases rapidly after a short delay.

☰ PARAMETER/MODE KEY

Used to enter parameter selections, access operating modes and index through menu items.

CAUTION

Changing parameters in these menus will most likely result in malfunction of the controller and possible damage to the source and/or controller.

4.1.3 Error Message Display

Problems in the Model 2500 are displayed as error messages on the °C display. Reported errors as A/D converter problems, and problem SENSOR conditions.

The problem messages and needed actions are as follows:

| DISPLAY | PROBLEM | ACTION |
|---------|-----------------------------|-----------------------|
| Err.H | Open sensor | Check sensor & wiring |
| Err.L | Shorted sensor | Check sensor & wiring |
| — — — — | Sensor over- or under-range | Check sensor & wiring |
| Err.O | A/D error | Return to EOI |
| Err.J | A/D error | Return to EOI |

4.2 BASIC OPERATING PROCEDURE IN LOCAL CONTROL

4.2.1 Power ON

Turn power ON. All LED segments and indicators are illuminated to verify proper operation. The display then identifies the temperature display units and the software revision number.

After the controller has been energized for 3 seconds, the °C display will indicate the temperature of the source and the SET display will indicate the current temperature set- point (the set-point will be the same as when the controller was last turned off).

4.2.2 Setting Temperature

The source temperature can be changed by pressing the RAISE ▲ or LOWER ▼ keys and watching the temperature value on the SET display. When either key is held down, the value changes rapidly after a short delay.

4.2.3 Cool Down Procedure

Set the blackbody temperature to 50°C. After it reaches 50°C, turn the power off.

CAUTION

If you turn the system off without first letting it cool down, high heat may damage the system. First let it cool to a lower temperature, then shut it off.

4.3 Remote Control (Option /R RS232 Interface)

The optional RS232 interface allows bi-directional data transfer via a three conductor cable consisting of signal ground, receive input and transmit output. It is recommended that the cable length be less than fifty feet between computer terminal and the controller. Note that multiple instruments cannot be connected to the same port.

The communication parameters are fixed at 9600 Baud, no parity, 8 data bits and 1 stop bit. No NULL MODEM needed.

RS232 connector pin-out (DB9S, female)

pin 2 ----- Tx (transmit line out of controller)
pin 3 ----- Rx (receive in to controller)
pin 5 ----- COM (signal common)

CAUTION

COM is a signal ground only. DO NOT connect to earth ground or damage to the controller may result.

4.3.1 Temperature Control

To set the temperature of the source via RS232 interface, send the temperature setpoint command in the proper format as described below. The setpoint command format:

```
[START CHAR] [ID] [TYPE] [PARAM] [DATA] [CHKSUM] [END  
CHAR]  
$      0101  W      09      xx.xxx      yy  
<CR>
```

There are two fields you need to fill in.

[DATA] is a six-character ASCII representation of a numeric value. All ASCII representations of numeric values must be done using the characters '0'-'9' and '.'.

Note:

1. Valid number should be between 0 - 1250. This range varies for each Source.
2. No space character allowed in this field.
3. Must be six characters critically (including '.'). If the number you want to set is less than 6 characters, add zeros before or after to make it 6-character long. 020.00 is valid.

[CHKSUM] is a two-character Message Code Numbering System, representing the sum of the ASCII values of all the characters (excluding the START CHAR, the END CHAR, and the CHKSUM themselves) in the command. The sum is computed using the following formula:

$$\text{CHKSUM} = \text{SUM}(\text{All Characters}) \% 256$$

% is the modulus operator.

In Message Code Numbering System, the most significant digit is represented with the numbers 0-9 and the letters A-Z and the least significant digit is represented with the numbers 0-9. The numbers 0-9 have the same values as their decimal counterparts and the letters A-Z have the values of 100 - 350 inclusively in increments of 10.

| CHKSUM value | = | Decimal value |
|--------------|---|---------------|
| "00" | | 0 + 0 = 0 |
| "99" | | 90 + 9 = 99 |
| "A0" | | 100 + 0 = 10 |
| "B8" | | 110 + 8 = 118 |
| "P5" | | 250 + 5 = 255 |

After sending a setpoint command, you will receive a response message from the controller. You must read the response message after sending a command.

The response message has the format:

[START CHAR] [ID] [TYPE] [PARAM] [ERROR] [CHKSUM] [END CHAR]

% 0101 W 09 x yy <CR>

[ERROR] field shows if the communication was successfully.
The data in this field would be:

0 means No error. Then, the number in [DATA] field will display on the Front Panel.

3 means Transmission failed: Parity error.

5 means Transmission failed: Bad message. Message cannot be understood.

6 means Transmission failed: Bad checksum. The checksum received did not match the checksum of the message.

A means Transmission failed: Bad data. Bad representation in the [DATA] field or data is out of range.

[CHKSUM] is the same as above. You can use it as the way to check if you have successfully received the message.

EXAMPLES:

| <u>Command</u> | <u>Temp. Setting</u> | <u>Response Message</u> | <u>Notes</u> |
|-------------------|---------------------------------|-------------------------|--------------|
| \$0101W0910.123G7 | 10.123°C | %0101W090H8 | No error |
| \$0101W09020.00G2 | 20°C | %0101W090H8 | No error |
| \$0101W09002000G5 | 2000°C | %0101W09AJ5 | Out of range |
| \$0101W0901000B5 | 1000°C | %0101W09 | Message. |
| [DATA] | field is less than 6 characters | | |

4.3.2 Temperature Readout

The current temperature as shown on the front panel display can be read via RS232 interface. A read command must first be sent to the controller. The command format:

[START CHAR] [ID] [TYPE] [PARAM] [CHKSUM] [END CHAR]

\$ 0101 R 05 C1 <CR>

[CHKSUM] field is the same with above. It is fixed as "C1" for this command. So, the message sent to the controller to read the temperature is:

\$0101R05C1<CR>

The Read response message has the format:

[START CHAR] [ID] [TYPE] [PARAM] [DATA] [CHKSUM] [END CHAR]

% 0101 R 05 xxx.xxx yy <CR>

where xxx.xxx is the current temperature and, as before, yy is the checksum.

EXAMPLES:

| <u>Response Message</u> | <u>Current Temp.</u> |
|-------------------------|----------------------|
| %0101R05016.304L3<CR> | 16.304°C |
| %0101R05020.030K4<CR> | 20.030°C |

4.3.3 Example Code

a) CHECKSUM calculation in C Language

```
// CHKSUM.C

// This short program is to show how to calculate the CHECKSUM in the setpoint
// command and how to append the CHECKSUM and Carriage Return character
// to // the original message.

#include "stdio.h"
#include "stdlib.h"
#include "string.h"
unsigned char LUT[26];
void main()
{
char RequestStr[15]="$0101W0910.123"; // For instance, Set Point = 10.123
char SendOut[17]; // SendOut contains original message, CHECKSUM and <CR>
int i;
unsigned int sum=0, rem;

for( i=0; i<=9; i++) // Prepare the Look Up Table {
LUT[i] = 0x30 + i; // for 0 – 9 }
for( i=10; i<26; i++) {
LUT[i] = 0x41 + i - 10; // for A – P }

// Calculate the sum of characters excluding the START CHAR
for(i=1; i<strlen(RequestStr); i++)
sum += RequestStr[i]; // which is $ and the CHKSUM themselves.

rem = sum % 256; // get modulus when sum is divided by 256
```

```
// append the two characters CHKSUM and <CR>( 0x0d ) to the original  
message  
sprintf(SendOut,"%s%c%c%c",RequestStr,LUT[rem/10],LUT[rem-  
rem/10*10],'\x0d');  
  
printf("%s",SendOut); }
```

b) Communications Software in LabVIEW

See Remote.exe and Remote.llb(LabVIEW source code) on accompanying CD.

SECTION 5

THEORY OF OPERATION

5.1 GENERAL

Referring to Interconnect drawing for Model 2500 series, temperature measurement and control functions are provided by A2. It also provides the user input for set point and display functions. K1 is a zero-crossing solid state relay (SSR) that controls the power to the source driven by a pulse-width modulated (PWM) signal from A2.

5.2 TEMPERATURE MEASUREMENT AND CONTROL (A2)

This module reads the source thermocouple connected at terminals 9 and 10 and converts the reading to temperature for display on the front panel. The temperature set point is changed with the front panel keys.

The module provides a PWM drive signal at terminals 1 and 2 to SSR K1 to control the power to the source. The line voltage is stepped down by transformer T1 (Model 2500D only) to provide low voltage drive to the source.

SECTION 6

MAINTENANCE

WARNING

These servicing instructions are for qualified service personnel only. To avoid electrical shock, do not perform any servicing other than that contained in the operating instructions unless you are qualified to do so.

6.1 PREVENTIVE MAINTENANCE

Under normal conditions, no preventive maintenance procedures are necessary for operation of the temperature controller.

6.2 TROUBLESHOOTING

In a complex microprocessor-based instrument such as this, a specific troubleshooting guide cannot be generated. Repair of the instrument will require a thorough understanding of the circuitry and careful application of good troubleshooting techniques. A review of the THEORY OF OPERATION in Section 5, would be helpful before attempting repair.

Frequently, the symptoms point to failure of one specific card, and the most expedient repair would be replacement of the entire card.

If repair is necessary, it is strongly recommended that the instrument be returned to Electro Optical Industries, Inc. where appropriate test equipment, trained service personnel and a complete stock of replacement parts can assure a satisfactory repair.

SECTION 7

CALIBRATION

7.1 INTRODUCTION

Calibration of the blackbody source and temperature controller is performed at Electro Optical Industries, Inc. If for some reason the accuracy of the source has changed, it is recommended that the source and controller be returned to Electro Optical Industries for checkout and re-calibration.

The recommended calibration cycle for EOI sources is as follows:

- (1) 50°C to 600°C Range, 5000 hours.
- (2) 80°C to 1000°C Range, 2000 hours.
- (3) 500 K to 1500 K Range, 1000 hours.

This cycle coincides with the expected shift in the control sensor in the source. A calibration procedure is provided in the following section.

7.2 RECOMMENDED EQUIPMENT FOR CALIBRATION

- (1) Calibrated and certified Pt/Pt + 10% Rh thermocouple whose calibration is known at three or more data points (similar to EOI Model 153S).
- (2) A laboratory stand to hold the thermocouple during calibration.
- (3) EOI Model 506 Potentiometric Microvoltmeter or equivalent.
- (4) An ice bath dewar flask such as Thermos bottle (for thermocouple cold reference probe).

7.3 USE OF CALIBRATION DATA SHEETS (Refer to Figure 7.1)

Definition of terms used on calibration sheets.

- (1) Measured mV: output of a thermocouple in millivolts.
- (2) Measured Temperature in °C: "Measured mV" converted to temperature by using a conversion chart for type S thermocouples. (See Table 7.1)

-
- (3) Error Temperature in °C: deviation from the conversion chart (for type S thermocouple, see Table 7.1 at the end of this section) for a specific thermocouple. When using a calibrated thermocouple, plot the error (error vs. temp) on graph paper by using the error points marked on the thermocouple tag. You can then find the error to any temperature by interpolation.
- (4) True Temperature in °C: actual source temperature found by combining the thermocouple "error temp." in °C with the temperature found on the conversion chart for type S thermocouple.



| | |
|-----------------|------------------------------|
| Blackbody Model | Potentiometer # |
| Controller # | Certified Thermocouple Model |
| Contract # | Customer |
| Serial # | Identification |
| Blackbody Res. | NIST |

Date _____

Signed _____

16

7.4 CALIBRATION SET-UP PROCEDURE (Manually from front panel) (Figure 7.2)

- 1) Set up a laboratory ring stand with a thermometer clamp.
- 2) Mount the calibrated thermocouple header in the thermometer clamp.
- 3) Adjust the height of clamped thermocouple until it is on center with blackbody cavity.

CAUTION

DO NOT insert or remove the thermocouple while cavity is hot. Maximum safe temperature is 100°C.

- 4) Push thermocouple into blackbody cavity until the thermocouple probe tip is in direct contact with the apex of the cavity.
- 5) Thermocouples with cold junction probe, insert into a crushed ice bath. Without cold junction probe, attach extension leads to an electronic 0°C reference junction.
- 6) Connect thermocouple leads to a potentiometric microvoltmeter.

7.5 CALIBRATION PROCEDURE

- 1) Turn temperature controller power ON and use the keypad to select the low end temperature of 50°C.
- 2) Allow sufficient time for stabilization at this temperature setting.
- 3) Record the measured millivolts (mV) of the calibrated thermocouple.
- 4) Convert into temperature using the NIST thermocouple reference table (see Table 7.1 at the end of this section).
- 5) To determine the true temperature of the source, plot an error curve from the calibration points given with the calibrated thermocouple.
 - a) If the error is negative, add to the source temperature.

- b) If the error is positive, subtract from the source temperature.
- 6) Next, Set the source temperature to 100° C and allow sufficient time for the blackbody to stabilize.
- 7) Record the measured millivolt output of the calibrated thermocouple. Continue this process as in steps 6 & 7 for every 100° until you have completed the blackbody temperature calibration.
- 8) After completion of calibration, set the temperature of the source to 50° C and allow the blackbody to cool. Pull the thermocouple out approximately 1½ inch so it's not touching the apex of the cavity.

CAUTION

Do not remove the calibrated thermocouple until the blackbody has cooled down to a Temperature of 100°C or below.

- 9) Convert into temperature using the NIST thermocouple reference table (see Table 7.1 at the end of this section).
- 10) To determine the true temperature of the source, plot an error curve from the calibration points given with the calibrated thermocouple.
 - a) If the error is negative, add to the source temperature.
 - b) If the error is positive, subtract from the source temperature.




Figure 7.2 Calibration Equipment Set up

See Appendix for Type S Thermocouple Table

SECTION 8

REPLACEABLE PARTS

8.1 ORDERING INFORMATION

To obtain replacement parts, order from the address on the front of the manual. Identify parts by EOI STOCK CODE, if this can be determined.

For non-listed parts, include the instrument Model and Serial numbers, the location of the part (Interconnect number or Component designator) and a description of the part (e.g., 10 μ F 35 V Tantalum Capacitor).

Because of the complexity of our instruments, or the microelectronics in them, some parts are not available and repairs should, or must, be done in our facility.

For very old equipment, it is helpful to know the original purchaser name, the purchase order number and the purchase date. We keep records of everything we sell. We can investigate these records to satisfy your requests. The more information we have, the easier it is to find what you need. Some parts for old equipment are no longer available. In these cases, we can try to find compatible parts to get your system back in service.

Here is a list of information that will help us find your parts:

- Model number;
- Serial number;
- Interconnect Number and revision letter;
- Component name and description;
- Any know specifications of component;
- Original purchaser name;
- Original purchase date;
- Purchase order number;
- Information on what the problem is, how often and under what circumstances it occurs.

LIMITED WARRANTY

Electro Optical Industries, Inc. (EOI) warrants this new Product against defects in materials or workmanship as follows:

For a period of Two (2) years from the date of shipment, EOI will repair or replace the Product, at EOI's option, without charge for parts or labor. After the period of Two (2) years their purchaser must pay all parts and labor charges.

The limited warranty is extended only to the original purchaser. It does not cover damages or failure caused by or attributed to Acts of God, abuse, misuse, improper or abnormal usage, faulty installation, improper maintenance, war, lightning, or other incidences of excessive voltage, or any repairs or tempering by other than an EOI authorized repair facility. It does not cover replacement of consumable parts, transportation costs, or damage in transit. *For the purposes of this warranty heater elements, emitting surface coatings and re-calibration are considered consumables.* This warranty will become void if any factory seals are removed or broken, or the serial number or model number identification has been wholly or partially removed or erased. Repair or replacement under the terms of this warranty does not extend the terms of this warranty. The specific components which are repaired or replaced are warranted for a period of six months. This warranty can not be modified by any agent of EOI, except in writing and signed by an officer of EOI.

Should this product prove defective in workmanship or material, the purchaser's sole remedies shall be such repair or replacement as provided by the terms of this warranty. Under no circumstances shall EOI be liable for any loss or damage, direct, consequential, or incidental, arising out of the use of or inability to use this Product.

To obtain warranty services, call or write EOI for a Return Material Authorization (RMA) number. Packaging must be suitable for the form of shipment used. Mark the RMS number clearly on the outside of the package. Include a copy of a description of the symptoms, contact name, organization, address, phone number and any special return shipping instructions. Then deliver or ship the product, postage or shipping costs prepaid, to EOI.

If Product is returned to EOI during the limited warranty period, but problem with the Product is determined to not be covered under the terms and conditions of this limited warranty the purchaser will be notified and given an estimate of the charges the

purchaser must pay to have the Product repaired, with all shipping charges billed to the purchaser. If the estimate is refused, the Product will be returned COD for the standard evaluation charge, if any, freight collect. If the Product is returned to EOI after the expiration of the warranty period, EOI's normal service policies shall apply and the purchaser will be invoiced for all shipping charges.

Purchaser may also have other rights, which may vary from state to state.

ANY IMPLIED WARRANTY OF MERCHANTABILITY, OR FITNESS FOR A PARTICULAR PURPOSE OR USE, SHALL BE LIMITED TO THE DURATION OF THE FOREGOING WRITTEN WARRANTY, OTHERWISE, THE FOREGOING WARRANTY IS THE PURCHASER'S SOLE AND EXCLUSIVE REMEDY AND IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. EOI SHALL NOT BE LIABLE FOR INCIDENTAL OR CONSEQUENTIAL DAMAGES OR LOSS OF ANTICIPATED BENEFITS OR PROFITS, RESULTING FROM THE PURCHASE OR USE OF THE PRODUCT ARISING FROM THE BREACH OF THE WARRANTY EVEN IF EOI KNEW OF THE LIKELIHOOD OF SUCH DAMAGES.

* * *

APPENDIX
TABLE 7

**Type S Thermocouples - Thermoelectric voltage as a function of temperature (°C)
Reference Tables N.I.S.T. Monograph 175 Revised to ITS-90.**

| Thermoelectric Voltage in millivolts for Type S Thermocouples | | | | | | | | | | | | |
|---|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|
| deg C | 0 | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | deg C |
| 0 | 0.000 | 0.005 | 0.011 | 0.016 | 0.022 | 0.027 | 0.033 | 0.038 | 0.044 | 0.050 | 0.055 | 0 |
| 10 | 0.055 | 0.061 | 0.067 | 0.072 | 0.078 | 0.084 | 0.090 | 0.095 | 0.101 | 0.107 | 0.113 | 10 |
| 20 | 0.113 | 0.119 | 0.125 | 0.131 | 0.137 | 0.143 | 0.149 | 0.155 | 0.161 | 0.167 | 0.173 | 20 |
| 30 | 0.173 | 0.179 | 0.185 | 0.191 | 0.197 | 0.204 | 0.210 | 0.216 | 0.222 | 0.229 | 0.235 | 30 |
| 40 | 0.235 | 0.241 | 0.248 | 0.254 | 0.260 | 0.267 | 0.273 | 0.280 | 0.286 | 0.292 | 0.299 | 40 |
| 50 | 0.299 | 0.305 | 0.312 | 0.319 | 0.325 | 0.332 | 0.338 | 0.345 | 0.352 | 0.358 | 0.365 | 50 |
| 60 | 0.365 | 0.372 | 0.378 | 0.385 | 0.392 | 0.399 | 0.405 | 0.412 | 0.419 | 0.426 | 0.433 | 60 |
| 70 | 0.433 | 0.440 | 0.446 | 0.453 | 0.460 | 0.467 | 0.474 | 0.481 | 0.488 | 0.495 | 0.502 | 70 |
| 80 | 0.502 | 0.509 | 0.516 | 0.523 | 0.530 | 0.538 | 0.545 | 0.552 | 0.559 | 0.566 | 0.573 | 80 |
| 90 | 0.573 | 0.580 | 0.588 | 0.595 | 0.602 | 0.609 | 0.617 | 0.624 | 0.631 | 0.639 | 0.646 | 90 |
| 100 | 0.646 | 0.653 | 0.661 | 0.668 | 0.675 | 0.683 | 0.690 | 0.698 | 0.705 | 0.713 | 0.720 | 100 |
| 110 | 0.720 | 0.727 | 0.735 | 0.743 | 0.750 | 0.758 | 0.765 | 0.773 | 0.780 | 0.788 | 0.795 | 110 |
| 120 | 0.795 | 0.803 | 0.811 | 0.818 | 0.826 | 0.834 | 0.841 | 0.849 | 0.857 | 0.865 | 0.872 | 120 |
| 130 | 0.872 | 0.880 | 0.888 | 0.896 | 0.903 | 0.911 | 0.919 | 0.927 | 0.935 | 0.942 | 0.950 | 130 |
| 140 | 0.950 | 0.958 | 0.966 | 0.974 | 0.982 | 0.990 | 0.998 | 1.006 | 1.013 | 1.021 | 1.029 | 140 |
| 150 | 1.029 | 1.037 | 1.045 | 1.053 | 1.061 | 1.069 | 1.077 | 1.085 | 1.094 | 1.102 | 1.110 | 150 |
| 160 | 1.110 | 1.118 | 1.126 | 1.134 | 1.142 | 1.150 | 1.158 | 1.167 | 1.175 | 1.183 | 1.191 | 160 |
| 170 | 1.191 | 1.199 | 1.207 | 1.216 | 1.224 | 1.232 | 1.240 | 1.249 | 1.257 | 1.265 | 1.273 | 170 |
| 180 | 1.273 | 1.282 | 1.290 | 1.298 | 1.307 | 1.315 | 1.323 | 1.332 | 1.340 | 1.348 | 1.357 | 180 |
| 190 | 1.357 | 1.365 | 1.373 | 1.382 | 1.390 | 1.399 | 1.407 | 1.415 | 1.424 | 1.432 | 1.441 | 190 |
| 200 | 1.441 | 1.449 | 1.458 | 1.466 | 1.475 | 1.483 | 1.492 | 1.500 | 1.509 | 1.517 | 1.526 | 200 |
| 210 | 1.526 | 1.534 | 1.543 | 1.551 | 1.560 | 1.569 | 1.577 | 1.586 | 1.594 | 1.603 | 1.612 | 210 |
| 220 | 1.612 | 1.620 | 1.629 | 1.638 | 1.646 | 1.655 | 1.663 | 1.672 | 1.681 | 1.690 | 1.698 | 220 |
| 230 | 1.698 | 1.707 | 1.716 | 1.724 | 1.733 | 1.742 | 1.751 | 1.759 | 1.768 | 1.777 | 1.786 | 230 |
| 240 | 1.786 | 1.794 | 1.803 | 1.812 | 1.821 | 1.829 | 1.838 | 1.847 | 1.856 | 1.865 | 1.874 | 240 |
| 250 | 1.874 | 1.882 | 1.891 | 1.900 | 1.909 | 1.918 | 1.927 | 1.936 | 1.944 | 1.953 | 1.962 | 250 |

| Thermoelectric Voltage in millivolts for Type S Thermocouples | | | | | | | | | | | | |
|---|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|
| deg C | 0 | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | deg C |
| 260 | 1.962 | 1.971 | 1.980 | 1.989 | 1.998 | 2.007 | 2.016 | 2.025 | 2.034 | 2.043 | 2.052 | 260 |
| 270 | 2.052 | 2.061 | 2.070 | 2.078 | 2.087 | 2.096 | 2.105 | 2.114 | 2.123 | 2.132 | 2.141 | 270 |
| 280 | 2.141 | 2.151 | 2.160 | 2.169 | 2.178 | 2.187 | 2.196 | 2.205 | 2.214 | 2.223 | 2.232 | 280 |
| 290 | 2.232 | 2.241 | 2.250 | 2.259 | 2.268 | 2.277 | 2.287 | 2.296 | 2.305 | 2.314 | 2.323 | 290 |
| 300 | 2.323 | 2.332 | 2.341 | 2.350 | 2.360 | 2.369 | 2.378 | 2.387 | 2.396 | 2.405 | 2.415 | 300 |
| 310 | 2.415 | 2.424 | 2.433 | 2.442 | 2.451 | 2.461 | 2.470 | 2.479 | 2.488 | 2.497 | 2.507 | 310 |
| 320 | 2.507 | 2.516 | 2.525 | 2.534 | 2.544 | 2.553 | 2.562 | 2.571 | 2.581 | 2.590 | 2.599 | 320 |
| 330 | 2.599 | 2.609 | 2.618 | 2.627 | 2.636 | 2.646 | 2.655 | 2.664 | 2.674 | 2.683 | 2.692 | 330 |
| 340 | 2.692 | 2.702 | 2.711 | 2.720 | 2.730 | 2.739 | 2.748 | 2.758 | 2.767 | 2.776 | 2.786 | 340 |
| 350 | 2.786 | 2.795 | 2.805 | 2.814 | 2.823 | 2.833 | 2.842 | 2.851 | 2.861 | 2.870 | 2.880 | 350 |
| 360 | 2.880 | 2.889 | 2.899 | 2.908 | 2.917 | 2.927 | 2.936 | 2.946 | 2.955 | 2.965 | 2.974 | 360 |
| 370 | 2.974 | 2.983 | 2.993 | 3.002 | 3.012 | 3.021 | 3.031 | 3.040 | 3.050 | 3.059 | 3.069 | 370 |
| 380 | 3.069 | 3.078 | 3.088 | 3.097 | 3.107 | 3.116 | 3.126 | 3.135 | 3.145 | 3.154 | 3.164 | 380 |
| 390 | 3.164 | 3.173 | 3.183 | 3.192 | 3.202 | 3.212 | 3.221 | 3.231 | 3.240 | 3.250 | 3.259 | 390 |
| 400 | 3.259 | 3.269 | 3.279 | 3.288 | 3.298 | 3.307 | 3.317 | 3.326 | 3.336 | 3.346 | 3.355 | 400 |
| 410 | 3.355 | 3.365 | 3.374 | 3.384 | 3.394 | 3.403 | 3.413 | 3.423 | 3.432 | 3.442 | 3.451 | 410 |
| 420 | 3.451 | 3.461 | 3.471 | 3.480 | 3.490 | 3.500 | 3.509 | 3.519 | 3.529 | 3.538 | 3.548 | 420 |
| 430 | 3.548 | 3.558 | 3.567 | 3.577 | 3.587 | 3.596 | 3.606 | 3.616 | 3.626 | 3.635 | 3.645 | 430 |
| 440 | 3.645 | 3.655 | 3.664 | 3.674 | 3.684 | 3.694 | 3.703 | 3.713 | 3.723 | 3.732 | 3.742 | 440 |
| 450 | 3.742 | 3.752 | 3.762 | 3.771 | 3.781 | 3.791 | 3.801 | 3.810 | 3.820 | 3.830 | 3.840 | 450 |
| 460 | 3.840 | 3.850 | 3.859 | 3.869 | 3.879 | 3.889 | 3.898 | 3.908 | 3.918 | 3.928 | 3.938 | 460 |
| 470 | 3.938 | 3.947 | 3.957 | 3.967 | 3.977 | 3.987 | 3.997 | 4.006 | 4.016 | 4.026 | 4.036 | 470 |
| 480 | 4.036 | 4.046 | 4.056 | 4.065 | 4.075 | 4.085 | 4.095 | 4.105 | 4.115 | 4.125 | 4.134 | 480 |
| 490 | 4.134 | 4.144 | 4.154 | 4.164 | 4.174 | 4.184 | 4.194 | 4.204 | 4.213 | 4.223 | 4.233 | 490 |
| 500 | 4.233 | 4.243 | 4.253 | 4.263 | 4.273 | 4.283 | 4.293 | 4.303 | 4.313 | 4.323 | 4.332 | 500 |
| 510 | 4.332 | 4.342 | 4.352 | 4.362 | 4.372 | 4.382 | 4.392 | 4.402 | 4.412 | 4.422 | 4.432 | 510 |
| 520 | 4.432 | 4.442 | 4.452 | 4.462 | 4.472 | 4.482 | 4.492 | 4.502 | 4.512 | 4.522 | 4.532 | 520 |
| 530 | 4.532 | 4.542 | 4.552 | 4.562 | 4.572 | 4.582 | 4.592 | 4.602 | 4.612 | 4.622 | 4.632 | 530 |
| 540 | 4.632 | 4.642 | 4.652 | 4.662 | 4.672 | 4.682 | 4.692 | 4.702 | 4.712 | 4.722 | 4.732 | 540 |
| 550 | 4.732 | 4.742 | 4.752 | 4.762 | 4.772 | 4.782 | 4.793 | 4.803 | 4.813 | 4.823 | 4.833 | 550 |
| 560 | 4.833 | 4.843 | 4.853 | 4.863 | 4.873 | 4.883 | 4.893 | 4.904 | 4.914 | 4.924 | 4.934 | 560 |

| Thermoelectric Voltage in millivolts for Type S Thermocouples | | | | | | | | | | | | |
|---|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|
| deg C | 0 | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | deg C |
| 570 | 4.934 | 4.944 | 4.954 | 4.964 | 4.974 | 4.984 | 4.995 | 5.005 | 5.015 | 5.025 | 5.035 | 570 |
| 580 | 5.035 | 5.045 | 5.055 | 5.066 | 5.076 | 5.086 | 5.096 | 5.106 | 5.116 | 5.127 | 5.137 | 580 |
| 590 | 5.137 | 5.147 | 5.157 | 5.167 | 5.178 | 5.188 | 5.198 | 5.208 | 5.218 | 5.228 | 5.239 | 590 |
| 600 | 5.239 | 5.249 | 5.259 | 5.269 | 5.280 | 5.290 | 5.300 | 5.310 | 5.320 | 5.331 | 5.341 | 600 |
| 610 | 5.341 | 5.351 | 5.361 | 5.372 | 5.382 | 5.392 | 5.402 | 5.413 | 5.423 | 5.433 | 5.443 | 610 |
| 620 | 5.443 | 5.454 | 5.464 | 5.474 | 5.485 | 5.495 | 5.505 | 5.515 | 5.526 | 5.536 | 5.546 | 620 |
| 630 | 5.546 | 5.557 | 5.567 | 5.577 | 5.588 | 5.598 | 5.608 | 5.618 | 5.629 | 5.639 | 5.649 | 630 |
| 640 | 5.649 | 5.660 | 5.670 | 5.680 | 5.691 | 5.701 | 5.712 | 5.722 | 5.732 | 5.743 | 5.753 | 640 |
| 650 | 5.753 | 5.763 | 5.774 | 5.784 | 5.794 | 5.805 | 5.815 | 5.826 | 5.836 | 5.846 | 5.857 | 650 |
| 660 | 5.857 | 5.867 | 5.878 | 5.888 | 5.898 | 5.909 | 5.919 | 5.930 | 5.940 | 5.950 | 5.961 | 660 |
| 670 | 5.961 | 5.971 | 5.982 | 5.992 | 6.003 | 6.013 | 6.024 | 6.034 | 6.044 | 6.055 | 6.065 | 670 |
| 680 | 6.065 | 6.076 | 6.086 | 6.097 | 6.107 | 6.118 | 6.128 | 6.139 | 6.149 | 6.160 | 6.170 | 680 |
| 690 | 6.170 | 6.181 | 6.191 | 6.202 | 6.212 | 6.223 | 6.233 | 6.244 | 6.254 | 6.265 | 6.275 | 690 |
| 700 | 6.275 | 6.286 | 6.296 | 6.307 | 6.317 | 6.328 | 6.338 | 6.349 | 6.360 | 6.370 | 6.381 | 700 |
| 710 | 6.381 | 6.391 | 6.402 | 6.412 | 6.423 | 6.434 | 6.444 | 6.455 | 6.465 | 6.476 | 6.486 | 710 |
| 720 | 6.486 | 6.497 | 6.508 | 6.518 | 6.529 | 6.539 | 6.550 | 6.561 | 6.571 | 6.582 | 6.593 | 720 |
| 730 | 6.593 | 6.603 | 6.614 | 6.624 | 6.635 | 6.646 | 6.656 | 6.667 | 6.678 | 6.688 | 6.699 | 730 |
| 740 | 6.699 | 6.710 | 6.720 | 6.731 | 6.742 | 6.752 | 6.763 | 6.774 | 6.784 | 6.795 | 6.806 | 740 |
| 750 | 6.806 | 6.817 | 6.827 | 6.838 | 6.849 | 6.859 | 6.870 | 6.881 | 6.892 | 6.902 | 6.913 | 750 |
| 760 | 6.913 | 6.924 | 6.934 | 6.945 | 6.956 | 6.967 | 6.977 | 6.988 | 6.999 | 7.010 | 7.020 | 760 |
| 770 | 7.020 | 7.031 | 7.042 | 7.053 | 7.064 | 7.074 | 7.085 | 7.096 | 7.107 | 7.117 | 7.128 | 770 |
| 780 | 7.128 | 7.139 | 7.150 | 7.161 | 7.172 | 7.182 | 7.193 | 7.204 | 7.215 | 7.226 | 7.236 | 780 |
| 790 | 7.236 | 7.247 | 7.258 | 7.269 | 7.280 | 7.291 | 7.302 | 7.312 | 7.323 | 7.334 | 7.345 | 790 |
| 800 | 7.345 | 7.356 | 7.367 | 7.378 | 7.388 | 7.399 | 7.410 | 7.421 | 7.432 | 7.443 | 7.454 | 800 |
| 810 | 7.454 | 7.465 | 7.476 | 7.487 | 7.497 | 7.508 | 7.519 | 7.530 | 7.541 | 7.552 | 7.563 | 810 |
| 820 | 7.563 | 7.574 | 7.585 | 7.596 | 7.607 | 7.618 | 7.629 | 7.640 | 7.651 | 7.662 | 7.673 | 820 |
| 830 | 7.673 | 7.684 | 7.695 | 7.706 | 7.717 | 7.728 | 7.739 | 7.750 | 7.761 | 7.772 | 7.783 | 830 |
| 840 | 7.783 | 7.794 | 7.805 | 7.816 | 7.827 | 7.838 | 7.849 | 7.860 | 7.871 | 7.882 | 7.893 | 840 |
| 850 | 7.893 | 7.904 | 7.915 | 7.926 | 7.937 | 7.948 | 7.959 | 7.970 | 7.981 | 7.992 | 8.003 | 850 |
| 860 | 8.003 | 8.014 | 8.026 | 8.037 | 8.048 | 8.059 | 8.070 | 8.081 | 8.092 | 8.103 | 8.114 | 860 |
| 870 | 8.114 | 8.125 | 8.137 | 8.148 | 8.159 | 8.170 | 8.181 | 8.192 | 8.203 | 8.214 | 8.226 | 870 |
| 880 | 8.226 | 8.237 | 8.248 | 8.259 | 8.270 | 8.281 | 8.293 | 8.304 | 8.315 | 8.326 | 8.337 | 880 |

| Thermoelectric Voltage in millivolts for Type S Thermocouples | | | | | | | | | | | | |
|---|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|-------|
| deg C | 0 | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | deg C |
| 890 | 8.337 | 8.348 | 8.360 | 8.371 | 8.382 | 8.393 | 8.404 | 8.416 | 8.427 | 8.438 | 8.449 | 890 |
| 900 | 8.449 | 8.460 | 8.472 | 8.483 | 8.494 | 8.505 | 8.517 | 8.528 | 8.539 | 8.550 | 8.562 | 900 |
| 910 | 8.562 | 8.573 | 8.584 | 8.595 | 8.607 | 8.618 | 8.629 | 8.640 | 8.652 | 8.663 | 8.674 | 910 |
| 920 | 8.674 | 8.685 | 8.697 | 8.708 | 8.719 | 8.731 | 8.742 | 8.753 | 8.765 | 8.776 | 8.787 | 920 |
| 930 | 8.787 | 8.798 | 8.810 | 8.821 | 8.832 | 8.844 | 8.855 | 8.866 | 8.878 | 8.889 | 8.900 | 930 |
| 940 | 8.900 | 8.912 | 8.923 | 8.935 | 8.946 | 8.957 | 8.969 | 8.980 | 8.991 | 9.003 | 9.014 | 940 |
| 950 | 9.014 | 9.025 | 9.037 | 9.048 | 9.060 | 9.071 | 9.082 | 9.094 | 9.105 | 9.117 | 9.128 | 950 |
| 960 | 9.128 | 9.139 | 9.151 | 9.162 | 9.174 | 9.185 | 9.197 | 9.208 | 9.219 | 9.231 | 9.242 | 960 |
| 970 | 9.242 | 9.254 | 9.265 | 9.277 | 9.288 | 9.300 | 9.311 | 9.323 | 9.334 | 9.345 | 9.357 | 970 |
| 980 | 9.357 | 9.368 | 9.380 | 9.391 | 9.403 | 9.414 | 9.426 | 9.437 | 9.449 | 9.460 | 9.472 | 980 |
| 990 | 9.472 | 9.483 | 9.495 | 9.506 | 9.518 | 9.529 | 9.541 | 9.552 | 9.564 | 9.576 | 9.587 | 990 |

| deg C | 0 | 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 | 10 | deg C |
|-------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|-------|
| 1000 | 9.587 | 9.599 | 9.610 | 9.622 | 9.633 | 9.645 | 9.656 | 9.668 | 9.680 | 9.691 | 9.703 | 1000 |
| 1010 | 9.703 | 9.714 | 9.726 | 9.737 | 9.749 | 9.761 | 9.772 | 9.784 | 9.795 | 9.807 | 9.819 | 1010 |
| 1020 | 9.819 | 9.830 | 9.842 | 9.853 | 9.865 | 9.877 | 9.888 | 9.900 | 9.911 | 9.923 | 9.935 | 1020 |
| 1030 | 9.935 | 9.946 | 9.958 | 9.970 | 9.981 | 9.993 | 10.005 | 10.016 | 10.028 | 10.040 | 10.051 | 1030 |
| 1040 | 10.051 | 10.063 | 10.075 | 10.086 | 10.098 | 10.110 | 10.121 | 10.133 | 10.145 | 10.156 | 10.168 | 1040 |
| 1050 | 10.168 | 10.180 | 10.191 | 10.203 | 10.215 | 10.227 | 10.238 | 10.250 | 10.262 | 10.273 | 10.285 | 1050 |
| 1060 | 10.285 | 10.297 | 10.309 | 10.320 | 10.332 | 10.344 | 10.356 | 10.367 | 10.379 | 10.391 | 10.403 | 1060 |
| 1070 | 10.403 | 10.414 | 10.426 | 10.438 | 10.450 | 10.461 | 10.473 | 10.485 | 10.497 | 10.509 | 10.520 | 1070 |
| 1080 | 10.520 | 10.532 | 10.544 | 10.556 | 10.567 | 10.579 | 10.591 | 10.603 | 10.615 | 10.626 | 10.638 | 1080 |
| 1090 | 10.638 | 10.650 | 10.662 | 10.674 | 10.686 | 10.697 | 10.709 | 10.721 | 10.733 | 10.745 | 10.757 | 1090 |
| 1100 | 10.757 | 10.768 | 10.780 | 10.792 | 10.804 | 10.816 | 10.828 | 10.839 | 10.851 | 10.863 | 10.875 | 1100 |
| 1110 | 10.875 | 10.887 | 10.899 | 10.911 | 10.922 | 10.934 | 10.946 | 10.958 | 10.970 | 10.982 | 10.994 | 1110 |
| 1120 | 10.994 | 11.006 | 11.017 | 11.029 | 11.041 | 11.053 | 11.065 | 11.077 | 11.089 | 11.101 | 11.113 | 1120 |
| 1130 | 11.113 | 11.125 | 11.136 | 11.148 | 11.160 | 11.172 | 11.184 | 11.196 | 11.208 | 11.220 | 11.232 | 1130 |
| 1140 | 11.232 | 11.244 | 11.256 | 11.268 | 11.280 | 11.291 | 11.303 | 11.315 | 11.327 | 11.339 | 11.351 | 1140 |
| 1150 | 11.351 | 11.363 | 11.375 | 11.387 | 11.399 | 11.411 | 11.423 | 11.435 | 11.447 | 11.459 | 11.471 | 1150 |
| 1160 | 11.471 | 11.483 | 11.495 | 11.507 | 11.519 | 11.531 | 11.542 | 11.554 | 11.566 | 11.578 | 11.590 | 1160 |
| 1170 | 11.590 | 11.602 | 11.614 | 11.626 | 11.638 | 11.650 | 11.662 | 11.674 | 11.686 | 11.698 | 11.710 | 1170 |
| 1180 | 11.710 | 11.722 | 11.734 | 11.746 | 11.758 | 11.770 | 11.782 | 11.794 | 11.806 | 11.818 | 11.830 | 1180 |
| 1200 | 11.951 | 11.963 | 11.975 | 11.987 | 11.999 | 12.011 | 12.023 | 12.035 | 12.047 | 12.059 | 12.071 | 1200 |
| 1210 | 12.071 | 12.083 | 12.095 | 12.107 | 12.119 | 12.131 | 12.143 | 12.155 | 12.167 | 12.179 | 12.191 | 1210 |
| 1220 | 12.191 | 12.203 | 12.216 | 12.228 | 12.240 | 12.252 | 12.264 | 12.276 | 12.288 | 12.300 | 12.312 | 1220 |
| 1230 | 12.312 | 12.324 | 12.336 | 12.348 | 12.360 | 12.372 | 12.384 | 12.397 | 12.409 | 12.421 | 12.433 | 1230 |
| 1240 | 12.433 | 12.445 | 12.457 | 12.469 | 12.481 | 12.493 | 12.505 | 12.517 | 12.529 | 12.542 | 12.554 | 1240 |
| 1250 | 12.554 | 12.566 | 12.578 | 12.590 | 12.602 | 12.614 | 12.626 | 12.638 | 12.650 | 12.662 | 12.675 | 1250 |
| 1260 | 12.675 | 12.687 | 12.699 | 12.711 | 12.723 | 12.735 | 12.747 | 12.759 | 12.771 | 12.783 | 12.796 | 1260 |
| 1270 | 12.796 | 12.808 | 12.820 | 12.832 | 12.844 | 12.856 | 12.868 | 12.880 | 12.892 | 12.905 | 12.917 | 1270 |
| 1280 | 12.917 | 12.929 | 12.941 | 12.953 | 12.965 | 12.977 | 12.989 | 13.001 | 13.014 | 13.026 | 13.038 | 1280 |
| 1290 | 13.038 | 13.050 | 13.062 | 13.074 | 13.086 | 13.098 | 13.111 | 13.123 | 13.135 | 13.147 | 13.159 | 1290 |
| 1310 | 13.280 | 13.292 | 13.305 | 13.317 | 13.329 | 13.341 | 13.353 | 13.365 | 13.377 | 13.390 | 13.402 | 1310 |
| 1320 | 13.402 | 13.414 | 13.426 | 13.438 | 13.450 | 13.462 | 13.474 | 13.487 | 13.499 | 13.511 | 13.523 | 1320 |

| | | | | | | | | | | | | |
|-------------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|-------------|
| 1330 | 13.523 | 13.535 | 13.547 | 13.559 | 13.572 | 13.584 | 13.596 | 13.608 | 13.620 | 13.632 | 13.644 | 1330 |
| 1340 | 13.644 | 13.657 | 13.669 | 13.681 | 13.693 | 13.705 | 13.717 | 13.729 | 13.742 | 13.754 | 13.766 | 1340 |
| 1350 | 13.766 | 13.778 | 13.790 | 13.802 | 13.814 | 13.826 | 13.839 | 13.851 | 13.863 | 13.875 | 13.887 | 1350 |
| 1360 | 13.887 | 13.899 | 13.911 | 13.924 | 13.936 | 13.948 | 13.960 | 13.972 | 13.984 | 13.996 | 14.009 | 1360 |
| 1370 | 14.009 | 14.021 | 14.033 | 14.045 | 14.057 | 14.069 | 14.081 | 14.094 | 14.106 | 14.118 | 14.130 | 1370 |
| 1380 | 14.130 | 14.142 | 14.154 | 14.166 | 14.178 | 14.191 | 14.203 | 14.215 | 14.227 | 14.239 | 14.251 | 1380 |
| 1390 | 14.251 | 14.263 | 14.276 | 14.288 | 14.300 | 14.312 | 14.324 | 14.336 | 14.348 | 14.360 | 14.373 | 1390 |
| | | | | | | | | | | | | |
| 1400 | 14.373 | 14.385 | 14.397 | 14.409 | 14.421 | 14.433 | 14.445 | 14.457 | 14.470 | 14.482 | 14.494 | 1400 |
| 1410 | 14.494 | 14.506 | 14.518 | 14.530 | 14.542 | 14.554 | 14.567 | 14.579 | 14.591 | 14.603 | 14.615 | 1410 |
| 1420 | 14.615 | 14.627 | 14.639 | 14.651 | 14.664 | 14.676 | 14.688 | 14.700 | 14.712 | 14.724 | 14.736 | 1420 |
| 1430 | 14.736 | 14.748 | 14.760 | 14.773 | 14.785 | 14.797 | 14.809 | 14.821 | 14.833 | 14.845 | 14.857 | 1430 |
| 1440 | 14.857 | 14.869 | 14.881 | 14.894 | 14.906 | 14.918 | 14.930 | 14.942 | 14.954 | 14.966 | 14.978 | 1440 |